



DAVID THOMAS

"an eye for excellence"

DTCL EQUIVALENT TO CIBA VISION RGP LENSES

Ciba Vision Gas permeable lenses

GAS PERMEABLE													
Minus BVP Range	PRODUCT NAME	DTCL PRODUCT NAME	BOZR mm	TD mm	H2O %	CT mm @-3.00	MATERIAL	TINT	dk	CYL	AXES	PRISM	Plus BVP Range
SPHERES & ASPHERES													
-30	POLYCON	FLUOROPERM 30, 60, 90, 151	7.00 to 8.60	9.00 9.50		0.1310.122	FLUOROCON 30, 60, 90, 151	Var.	12 to 151				30
-20	ASPHERICON	FLUOROPERM 30, 60, 90, 151	7.00 to 8.60	9.30 9.80		0.12	FLUOROCON 30, 60, 90, 151	Var.	12 to 95				20
-20	ASTROCON	FLUOROPERM 30, 60, 90, 151	7.00 to 8.60	9.60 10.20		0.135	FLUOROCON 30, 60, 90, 151, & HDS	Var.	12 to 95				20
-30	K DESIGN	FLUOROPERM 30, 60, 90, 151	7.00 to 8.60	9.00 9.50		0.17 0.15	FLUOROCON 30, 60, 90, 151	Var.	12 to 151				30
-30	KE DESIGN	FLUOROPERM 30, 60, 90, 151	7.00 to 8.60	9.20 9.80		0.15 0.16	FLUOROCON 30, 60, 90, 151	Var.	12 to 151				30
-30	PRACTITIONER SPECIFIED	FLUOROPERM 30, 60, 90, 151	to Rx	to Rx		to Rx	FLUOROCON 30, 60, 90, 151	Var.	12 to 151				30
-25	PERSECON® E	FLUOROLENS 25	6.90 to 8.60 7.00 to 8.70 7.10 to 8.80 7.20 to 8.90	8.80 9.30 9.80 10.30		0.17	FLUOROLENS 25	Var.	8				25
-25	PERSECON® 92E	FLUOROPERM 90	7.20 to 8.50 7.30 to 8.50 7.30 to 8.50	8.80 9.30 9.80 10.30		0.17	FLUOROCON 90	Var.	92				25
-25	PERSECON® CE	FLUOROPERM 60	7.20 to 8.40 7.30 to 8.50 7.30 to 8.50	8.80 9.30 9.80 10.30		0.17	FLUOROCON 60	Var.	54				25
-30	AQUILA	PARAGON HDS 100	7.40 TO 8.60	9.80 & 10.30		0.17	PARAGON HDS 100	Var.	100				30
-30	CONTINUOUS CURVE (PMMA)	PMMA	7.00 to 8.60	9.00 9.50		0.17 0.20	PMMA	(912) Grey	0				30
KERATOCONUS													
-30	SOPER KERATOCONUS	FLUOROPERM 30, 60	to Rx	to Rx		to Rx	FLUOROCON 30, 60, 90, 151	Var.	12 to 65				30
-30	WOODWARD KERATOCONUS	FLUOROPERM 30, 60	to Rx	to Rx		to Rx	FLUOROCON 30, 60, 90, 151	Var.	12 to 65				30
-25	PERSECON E KERATOCONUS	FLUOROLENS 25	6.00 to 7.70	9.30 9.80		0.203	FLUOROCON 30, 60, 90, 151	Clear	8				25
TORIC													
-30	FRONT SURFACE TORIC	FLUOROPERM 30, 60	7.00 to 8.60 (0.05)	8.00 to 10.00		to Rx	Various	Var.	12 to 65 Var.	up to 3.00D	1° ATC	up to 3.00Δ	30
-30	BACK SURFACE TORIC	FLUOROPERM 30, 60	7.00 to 8.60 (0.05)	8.00 to 10.00		to Rx	Various	Var.	12 to 65 Var.	up to 5.00D			30
-30	BI-TORIC	FLUOROPERM 30, 60	7.00 to 8.60 (0.05)	8.00 to 10.00		to Rx	Various	Var.	12 to 65 Var.	up to 5.00D			30
-30	TORIC PERIPHERY	FLUOROPERM 30, 60	to Rx	to Rx		to Rx	Various	Var.	12 to 65 Var.				30
MULTIFOCAL													
-30	ASTROCON	OBLONG GENIUS	6.50 TO 9.00	9.10 9.60 10.10		0.17	FLUOROCON 30, 60, 90, 151 & PARAGON HDS	Var.	12 to 151	ADD 1.00 TO 3.50D			30

PLEASE NOTE THAT THE FLUROCON RANGE OF MATERIALS SUPPLIED BY CIBA VISION ARE MADE FROM A RANGE OF MATERIALS CALLED FLUROPERM MANUFACTURED BY PARAGON VISION SCIENCES IN THE USA. UNLIKE SOME MANUFACTURERS DAVID THOMAS USES THIS EXACT SAME MATERIAL FROM PARAGON VISION SCIENCES.

PLEASE ALSO NOTE THAT DAVID THOMAS HAVE BEEN MANUFACTURING AND SUPPLYING THE ASTROCON MULTIFOCAL LENSES SINCE CIBA DIVESTED THIS BUSINESS APPROXIMATELY TWO YEARS AGO.

IF YOU REQUIRE AN ALTERNATIVE MATERIAL WE CAN SUPPLY THIS. THE CIBA RGP DESIGNS ARE NOT MATERIAL DEPENDENT.